
# I Etilometro ad alta precisione

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
# EN High-accuracy breathalyzer


- 1. LCD Display
- 2. Power button
- 3. Units of measurement
- 4. Battery compartment
- 5. Exhale pipe
- 6. Fixing the hand strap

## FEATURES

- Digital breathalyzer with audible alert.
- Sensor: advanced oxide alcohol semiconductor sensor.
- Range: 0.000-0.199% BAC (0.00-1.990‰ BAC o 0.000-0.995 mg/l BrAc)
- Accuracy: 0.005% BAC (0,1g/l).
- Backlight of the Led display in orange color.
- Quick response and resume.
- Auto power off: 15 seconds approx.
- Uses 2 x "AAA" alkaline batteries (not included).
- Operating time: warm up: 10-15 s, response <5 sec
- Dimensions: 116 x 45 x 15mm.
- Weight: 35gr.
- Environmental temperature: 10°C - 40°C.
- Environmental humidity: 40% - 90%.

## BATTERY INSTALLATION & REPLACEMENT


- 1 Open the battery compartment lid on the back of the unit.
- 2 Insert 2 "AAA" (UM-4) high-quality alkaline batteries (not included) with correct polarity as indicated.
- 3 Close the battery compartment lid.
- 4 When the LCD Display starts to dim, please replace old batteries.

## HOW TO USE THE BREATHALYZER

### Warm up

- 1 Press and hold the power button for 2 seconds, the LCD will turn on with one beeping sound.


- 2 The "WAIT" symbol is displayed, and the 10 to 0 countdown is started that indicate the tester is in warm-up mode.


**Note:** The sensor absorbs many impurities and gases present in the air when not in use. Therefore, it is advisable to repeat the warm-up phase several times before carrying out the test.

### TEST

- 1 The BLOW symbol is displayed on the screen, followed by a beep, indicating that the warm-up phase is finished and the alcohol tester is ready for use. The countdown from 10 to 0 will start, this means that the tester is in the testing mode.

- 2 Bring the tester closer and blow into the breath inhaler for about 3-5 seconds.

- 3 After a second beep it will be possible to read the result on the screen.

- 4 The alcohol tester will emit an acoustic signal if the alcohol concentration is equal or higher than 0,050% BAC (0,50‰ BAC o 0,250 mg/l BrAc).




- 5 If the tester does not detect any breath during the test phase (Blow), the value "0,000" will appear on the screen. To obtain a more reliable result, it is advisable to repeat the test 2 or 3 times.

- 4 The alcohol tester will emit an acoustic signal if the alcohol concentration is equal or higher than 0,050% BAC (0,50‰ BAC o 0,250 mg/l BrAc).

- 5 If the tester does not detect any breath during the test phase (Blow), the value "0,000" will appear on the screen. To obtain a more reliable result, it is advisable to repeat the test 2 or 3 times.


### UNITS OF MEASUREMENT

Press the unit button on the back of tester and select the measuring unit for the result: % BAC, ‰ BAC, mg/l BrAc.

## AUTO POWER OFF

The result is displayed on the screen for about 10 seconds after which the OFF symbol will appear and the alcohol tester will turn off.


## LOW BATTERY INDICATION

If the "LO" symbol flashes on the display, it indicates the batteries are almost empty and the tester will automatically turn off, please replace the old batteries.


## ERROR MESSAGE

If the "Err" symbol flashes on the display, it indicates a sensor error, and the tester will automatically turn off, please press again the power button and retry.

## ADJUSTING THE HAND STRAP

- 1 Insert the smaller loop into the hole in the housing.
- 2 Thread the longer loop through the eyelet and pull through.

## BACKLIGHT

When the tester is turned on, the LCD backlight will light until the tester auto powers off.

## PRECAUTIONS

After drinking, it is recommended to wait at least 20 minutes, as it takes approx. 20 minutes for alcohol to be absorbed into blood from the digestive organs, and residual alcohol remaining in the mouth takes this long to dissipate. For accurate result, do not repeat next testing within 3 minutes after last time. Do not smoke and eat for 20 minutes prior testing. Do not use chemicals such as benzene, thinner or petrol to clean the tester directly. Do not pour liquid into the breath pipe. Do not use any excitant substances near the tester, such as paints, insecticides, alcohols. Do not place the tester in a closed environment with contaminated air. Modifying or tampering with the unit's internal components may cause a malfunction. The unit is designed to be used as an aid to promote reasonable driving. It should not be used in anticipation of the use or operation of vehicles, machinery, or equipment. The unit makes no references whatsoever to blood alcohol content or other such measures of intoxication promulgated by government, medical or other agencies or bodies. The result of this alcohol gives BAC for reference only, manufacturer does not take any legal responsibility. Each individual's body responds differently to alcohol. A person may become legally impaired to drive at a point below the legal blood alcohol limit.

## DISCLAIMER

This tester indicates BAC (BrAc) for reference only. The manufacture will not take any legal responsibility for this.